

第 1 步. 使用 relion 估算 micrographs 的 CTF 参数, 生成 micrograph\_ctf.star

Step 1. Estimate CTF parameters for micrographs by RELION and generate *micrograph\_ctf.star*.

第 2 步. 将 star 文件转换为 GisSPA 可读的 lst 文件

Step 2. Convert *.star* files to GisSPA-readable *.lst* files

```
relion2lst_dont_write_mics.py CtfFind/job005/micrographs_ctf.star --lst micrographs_ctf.lst
```

```
#LST
0 MotionCorr/job003/micrographs/Yt_105kx_1202_00002.mrc defocus=0.8401005371 dfdiff=0.0128565918 dfang=145.699879
0 MotionCorr/job003/micrographs/Yt_105kx_1202_00003.mrc defocus=0.8660711426 dfdiff=0.0188164551 dfang=171.893494
0 MotionCorr/job003/micrographs/Yt_105kx_1202_00004.mrc defocus=0.86874335935 dfdiff=0.01908408205 dfang=170.796783
0 MotionCorr/job003/micrographs/Yt_105kx_1202_00005.mrc defocus=0.88477197265 dfdiff=0.02731279295 dfang=166.863052
0 MotionCorr/job003/micrographs/Yt_105kx_1202_00006.mrc defocus=0.84074106445 dfdiff=0.01255209965 dfang=174.929901
0 MotionCorr/job003/micrographs/Yt_105kx_1202_00007.mrc defocus=0.84785268555 dfdiff=0.01629379885 dfang=157.232521
0 MotionCorr/job003/micrographs/Yt_105kx_1202_00008.mrc defocus=0.87532441405 dfdiff=0.02016367185 dfang=156.715553
0 MotionCorr/job003/micrographs/Yt_105kx_1202_00009.mrc defocus=0.90362680665 dfdiff=0.02053491215 dfang=157.642876
0 MotionCorr/job003/micrographs/Yt_105kx_1202_00010.mrc defocus=0.84755024415 dfdiff=0.01133764645 dfang=140.132988
```

第 3 步. 为了节约计算时间, 可以先对图进行压缩, 根据像素大小可以选择 bin2

或者 bin4

Step 3. To improve computation efficiency, you can first compress the micrograph, and depending on the pixel size you can choose bin2 or bin4

```
cd MotionCorr/job003/micrographs/
```

```
mkdir ms2
```

```
foreach file (*.mrc)
e2proc2d.py $file ms2/$file --meanshrink=2
end
```

or

```
foreach file (*.mrc)
relion_image_handler --i $file --o ms2/$file --angpix 1.36 --rescale_angpix 2.72
end
```

替换 micrograph\_ctf.lst 的内容:

Replace the path of micrograph:

```
sed -i 's/micrographs/ms2\micrographs/g' micrographs_ctf.lst
```

```
#LST
0 MotionCorr/job003/ms2/micrographs/Yt_105kx_1202_00002.mrc defocus=0.8401005371 dfdiff=0.0128565918 dfang=145.699879
0 MotionCorr/job003/ms2/micrographs/Yt_105kx_1202_00003.mrc defocus=0.8660711426 dfdiff=0.0188164551 dfang=171.893494
0 MotionCorr/job003/ms2/micrographs/Yt_105kx_1202_00004.mrc defocus=0.86874335935 dfdiff=0.01908408205 dfang=170.796783
0 MotionCorr/job003/ms2/micrographs/Yt_105kx_1202_00005.mrc defocus=0.88477197265 dfdiff=0.02731279295 dfang=166.863052
0 MotionCorr/job003/ms2/micrographs/Yt_105kx_1202_00006.mrc defocus=0.84074106445 dfdiff=0.01255209965 dfang=174.929901
0 MotionCorr/job003/ms2/micrographs/Yt_105kx_1202_00007.mrc defocus=0.84785268555 dfdiff=0.01629379885 dfang=157.232521
0 MotionCorr/job003/ms2/micrographs/Yt_105kx_1202_00008.mrc defocus=0.87532441405 dfdiff=0.02016367185 dfang=156.715553
0 MotionCorr/job003/ms2/micrographs/Yt_105kx_1202_00009.mrc defocus=0.90362680665 dfdiff=0.02053491215 dfang=157.642876
0 MotionCorr/job003/ms2/micrographs/Yt_105kx_1202_00010.mrc defocus=0.84755024415 dfdiff=0.01133764645 dfang=140.132988
```

第 4 步. 将三维模型 rescale 成与 micrograph\_ctf.lst 相同的像素大小

**Step 4.** Rescale the 3D template to the same pixel size to the binned pixel size

```
relion_image_handler --i 3D_template.mrc --o 3D_template_rescaled.mrc --angpix ?? --rescale_angpix 2.72
```

**第 5 步.** 生成二维投影

**Step 5.** Generate 2D projections.

```
project3d 3D_template_rescaled.mrc out=2D_projections.hdf list=euler_C1_step5.lst
```

project3d 是 EMAN 软件包的一个用于创建二维投影的程序，下载 EMAN 即可使用 (<https://blake.bcm.edu/emanwiki/EMAN1>)。用于投影的欧拉角文件可在测试数据中下载。

*project3d* is a program for creating 2D projections from the EMAN package, EMAN software package can be downloaded from <https://blake.bcm.edu/emanwiki/EMAN1>. The Euler file can be downloaded from the test Data dir.

**第 6 步.** 写 config 文件。config 文件可在测试数据中下载，里面包含一些固定参数和可调参数

**Step 6.** Write config file, including some optional parameters.

```
input          = micrographs_ctf.lst
template       = 2D_projections.hdf
eulerfile      = euler_C1_step5.lst
angpix         = 2.72
phistep        = 3
kk             = 3
energy         = 300
cs             = 2.7
Highres        = 8
Lowres         = 400
diameter       = 160

norm_type      = 1
invert         = 1
threshold      = 6.2
output         = Output/test.lst
first          = 0
last           = 1
window_size    = 512
GPU_ID         = 0
phase_flip     = 1
overlap        = 160
```

input = input micrograph lstfile with ctf information

template = input 2D projections templates in .hdf format

eulerfile = euler file with euler values

angpix = input pixel size in angstroms

phistep = inplane rotation sampling  
 kk = overlapping density parameter, default is 3.  
 energy = accelerating voltage in kV.  
 cs = spherical aberration in um.  
 Highres = high resolution cut  
 Lowres = low resolution cut  
 diameter = target diameter in pixels Optional Parameters.  
 norm\_type = CCG norm optional parameter, 1 for global normalization, 0 for per-window normalization  
 invert = whether to invert contrast in micrograph, 1 for yes, 0 for no.  
 threshold = cc threshold value, only output score beyond this value  
 output = output lstfile filename  
 first = the first image id to process.  
 last = the last image id to process.  
 window\_size = the window size which is splitted from raw IMG, which must satisfy  $\text{window\_size} \% 32 == 0$ .  
 GPU\_ID = ID of GPU device, only support single ID.  
 phase\_flip = Whether do filtering on images, operation(1) or not(0), in case of input being filtered already).  
 overlap = size of overlap between different windows.

## 第 7 步. 运行程序进行颗粒探测

**Step 7.** Run program and do target detection.

```
~/GisSPA-main/main config
```

## 第 8 步. 合并重复探测的颗粒

**Step 8.** Merge duplicated particles

```
remove_repeat_particles_from_list.py Output/test.lst 1 4 8 Output/test_merge.lst
```

*remove\_repeat\_particles\_from\_list.py*

*<input1> <number of images> <center> <euler\_thres> <output>*

*<input1>* lst filename generated from target detection

*<number of images>* the number of micrographs used for detection in this input lst file

*<center>* the threshold of coordinate (x,y) range used to merge detected targets

*<euler\_thres>* the threshold of Euler angle (x,y) range used to merge detected targets

## 第 9 步. 将合并后的颗粒文件格式转换为 relion 可读的 star 文件

**Step 9.** Convert the merged lst file to the RELION-readable star file

```
convert_my_lst_to_relion_mic.py Output/test_merge.lst 2 1.36 test_merged.star MotionCorr/job003/micrographs/
```

*convert\_my\_lst\_to\_relion\_mic.py*

<input lstfile> <scale factor> <pixel size> <out star file> <micrographs dir>

<input lstfile> merged lst file

<scale factor> binning factor when target detection, 2 for bin2


<pixel size> the un-binned pixel size of original micrograph

<out star file> filename of output star file

<micrographs dir> the directory of original micrographs.


**第 10 步. 提取颗粒。使用 relion 的 re-extract 功能从 micrograph 中提取颗粒**



**Step 10. Particle extraction by RELION.**



**第 11 步. 三维分类**

**Step 11. 3D classification**





后续的操作和分析都与单颗粒的优化类似，但是必须是 local refinement，一定不能跑全局 refinement !!!

The subsequent operations and analysis are similar to the single particle analysis, but only local refinement is allowed, never try global refinement !!!

---

---

如果需要消除模型偏差，可使用 *jalign* 计算颗粒与重构之间的相位残差，计算的频率范围为：低频为 8 埃，高频由重构的分辨率决定。*jalign* 是 JSPR 软件包的一个程序，JSPR 下载链接为： <https://jiang.bio.purdue.edu/jspr/>

In case of model bias, the phase residuals between the particles and the reconstruction can be calculated using *jalign* for a frequency range of 8Å for low frequency cut and high frequency cut is determined by the resolution of the reconstruction. *jalign* is a program from the JSPR package, the JSPR can be downloaded at: <https://jiang.bio.purdue.edu/jspr/>

(1) 首先将 run\_data.star 转换为 lst 格式的文件，注意需要将 subset 1 和 subset 2 分开

(1) First convert run\_data.star to a file in lst format, note that subset 1 and subset 2 need to be separated

```
relion_star_handler --i run_data.star --o run_data_subset1.star --select rlnRandomSubset --minval 0.9 --maxval 1.1
relion_star_handler --i run_data.star --o run_data_subset2.star --select rlnRandomSubset --minval 1.9 --maxval 2.1
relion2lst_dont_write_particles.py run_data_subset1.star --lst run_data_subset1.lst --boxsize ???
relion2lst_dont_write_particles.py run_data_subset2.star --lst run_data_subset2.lst --boxsize ???
```

(2) 计算相位残差

(2) calculate phase residuals

```
$PATH/jalign run_data_subset2.lst run_class001.mrc out.lst --first $firstparticleID --last $lastparticleID
--apix $pixsize --aligner refineMicrographDefocus:batchsize=1:defocusrange=0.001:maskradius=$targetradius
:masksoft=-4:precision=0.001:stepdefocus=0.001:useOrigDefocus=1:verbose=1 --cmp phaseResRange:
ampweight=0:ctfampweight=1:lowRes=8:highRes=$resolution:abs=0 --batchsize -1 --force 1
```

(3) 选择合适的阈值，保留相关颗粒后重新转回 relion 格式

(3) Select particles with score above one threshold and re-convert it to star format

```
delete_bad_score.py $jspr_scored_lstfile $score_threshold $output_filename
jspr_refine_2_relion_class3d.py $selected_particles_lst $box_size $pixel_size $output_star_file 0/1 0/1/2
```